

Recombinant Murine IL-10

CATALOG #:	4156-10	10 µg
	4156-1000	1 mg
LOT #:	_____	
SOURCE:	<i>E. coli</i>	
PURITY:	>98% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of IL-10.	
MOL. WEIGHT:	18.7 kDa	
FORMULATION:	Lyophilized from 3 mM Tris, pH 8.0	

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in 3 mM Tris, pH 8.0 to a concentration of 0.1-1.0 µg/µl. The solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

STORAGE CONDITIONS:

The lyophilized IL-10 is best-stored desiccated below 0°C. Reconstituted IL-10 should be stored at working aliquots at -20°C.

DESCRIPTION:

IL-10 is an immunosuppressive cytokine produced by a variety of mammalian cell types including macrophages, monocytes, T cells, B cells and keratinocytes. IL-10 inhibits the expression of proinflammatory cytokines such as IL-1 and TNF- α . Like IL-4, IL-10 enhances humoral immune responses and attenuates cell-mediated immune reactions. Human IL-10 active on murine cells, but murine IL-10 is inactive on human cells. Recombinant murine IL-10 is an 18.7 kDa protein of 161 amino acid residues.

BIOLOGICAL ACTIVITY:

The ED₅₀ determined by the dose-dependent co-stimulation (with murine IL-4) of MC/9 cells is ≤ 2 ng/ml, corresponding to a specific activity of $\geq 5 \times 10^5$ units/mg.

USAGE: For Research Use Only. Not to be used in humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies