

Recombinant Human Interleukin-12 (IL-12)

CATALOG #:	4161-10	10 µg
	4161-1000	1 mg
LOT #:	_____	
SOURCE:	CHO cells (Chinese Hamster Ovarian Cells)	
PURITY:	>98% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of IL-12.	
MOL. WEIGHT:	75.0 kDa	
FORMULATION:	Sterile filtered and lyophilized in 0.5X PBS, pH 7.2	

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in 1X PBS (pH 7.2-7.4) to a concentration of 0.2 mg/ml. The solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

STORAGE CONDITIONS:

The lyophilized human IL-12 is best stored desiccated below 0°C. Reconstituted IL-12 should be stored at working aliquots at -20°C.

DESCRIPTION:

IL-12 is a potent regulator of cell mediated immune responses and induces IFN- γ production by NK and T cells. It is produced by activated monocytes/macrophage cells, B-lymphocytes, and connective tissue type mast cells. Among its biological activities, IL-12 promotes the growth and activity of activated NK, CD4+, and CD8+ cells and induce the development of IFN- γ producing Th1 cells. Recombinant human IL-12 is a 75 kDa heterodimeric glycoprotein consisting of disulfide-linked 35 kDa (p35) and 40 kDa (p40) subunits (503 total amino acid residues).

BIOLOGICAL ACTIVITY:

The ED₅₀ as determined by the induction of IFN- γ production from murine splenocytes co-stimulated with IL-12 is ≤ 0.1 ng/ml, corresponding to a specific activity of $\geq 1 \times 10^7$ units/mg.

USAGE: For Research Use Only. Not to be used in Humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies