

Recombinant Human Interleukin-13 (IL-13)

CATALOG #: 4164-10 10 µg
4164-1000 1 mg

LOT #: _____

SOURCE: *E. coli*

PURITY: >98% by SDS-PAGE and HPLC analysis
Endotoxin level is <0.1 ng per µg of IL-13.

MOL. WEIGHT: 12.5 kDa

FORMULATION: Sterile filtered and lyophilized from 20 mM HCl

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in 5-10 mM HCl to a concentration of 0.1 µg/µl. The solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

STORAGE CONDITIONS:

The lyophilized IL-13 should be stored desiccated below 0°C. Reconstituted IL-13 should be stored at working aliquots at -20°C.

DESCRIPTION:

IL-13 is an immunoregulatory protein produced by activated T lymphocytes. The biological activities of IL-13 include the stimulation of B cell proliferation, immunoglobulin production, and induction of differentiation of B cells and monocytic lineage. Recombinant human IL-13 is a 12.5 kDa protein containing 114 amino acid residues.

BIOLOGICAL ACTIVITY:

The ED₅₀ determined by the dose-dependent stimulation of the proliferation of TF-1 cells was ≤ 1.0 ng/ml, corresponding to a specific activity of ≥ 1 x 10⁶ units/mg.

USAGE: *For Research Use Only. Not to be used in human Subjects!*

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies