

Recombinant Murine IL-16

CATALOG #: 4174-10 10 µg
4174-1000 1 mg

LOT #: _____

SOURCE: *E. coli*

PURITY: >98% by SDS-PAGE and HPLC analysis
Endotoxin level is <0.1 ng per µg of IL-16.

MOL. WEIGHT: 13.2 kDa

FORMULATION: Lyophilized from 5 mM Sodium Phosphate, pH 7.6

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1 mg/ml. The solution can then be diluted into other aqueous buffers and store at 4°C for 1 week or – 20°C for future use.

STORAGE CONDITIONS:

The lyophilized IL-16 is best-stored desiccated below 0°C. Reconstituted IL-16 should be stored at working aliquots at -20°C.

DESCRIPTION:

Interleukin 16 (IL-16), also known as lymphocyte chemoattractant factor (LCF), stimulates a migratory response in CD4+ T cells and CD4+ monocytes and eosinophils. In addition, human IL-16 induces expression of the IL-2 receptor by T cells. Murine IL-16 is a 13.2 kDa protein containing 127 amino acid residues.

USAGE: For Research Purpose Only. Not to be used in humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies